ENGLISH

INION GTR™ BIODEGRADABLE MEMBRANE INSTRUCTIONS FOR USE

DESCRIPTION

The **INION GTR™ BIODEGRADABLE MEMBRANE** is intended to be used as a barrier membrane in dental guided tissue regeneration (GTR) and guided bone regeneration (GBR) procedures.

The **INION GTR™ BIODEGRADABLE SYSTEM** consists of membranes and separately available tacks made of degradable copolymers composed of L-lactic acid, D-lactic acid, glycolic acid and trimethylene carbonate. These polymers have a long history of safe medical use and they degrade in vivo by hydrolysis into alpha-hydroxy acids that are metabolised by the body. The membranes are fastened in position by use of separately available **INION GTR™ BIODEGRADABLE TACKS** or commercially available resorbable sutures with a non-cutting needle. The membrane provides a barrier over 8-12 weeks in vivo. Bioresorption takes place within one to two years. NMP (N-methyl-2-pyrrolidone) plasticizer temporarily softens the membrane and is used to facilitate the handling and shaping of the membrane. NMP is a pharmaceutical excipient that is used in several approved FDA parenteral drug products. NMP has a low degree of toxicity, is non-carcinogenic, non-mutagenic and undergoes rapid biodegradation.

The **INION GTR™ BIODEGRADABLE MEMBRANE** is designed to be used with customized instrumentation, e.g. the **INION GTR™** drill, tacks, tack applicators and shaping templates.

The **INION GTR™ BIODEGRADABLE SYSTEM** implants are sterile and non-collagenous.

INDICATIONS

The **INION GTR™ BIODEGRADABLE MEMBRANE** is indicated for the surgical treatment of periodontal defects to aid in the regeneration and integration of tissue components in guided tissue regeneration procedures (e.g. class II furcation defects, intrabony defects and recession type defects), for pre-implant and peri-implant surgery and for covering bone defects and empty sockets.

CONTRAINDICATIONS

The **INION GTR™ BIODEGRADABLE MEMBRANE** should not be used in patients with active or potential infections, patient conditions including limited blood supply, when patient cooperation cannot be guaranteed (e.g., alcoholism, drug abuse), and in load bearing indications unless used in conjunction with traditional rigid fixation. There are currently no known additional contraindications to the use of the **INION GTR™ BIODEGRADABLE MEMBRANE**.

INFORMATION FOR USE

Surgical considerations and reminders

- Prophylactic preoperative antibiotic treatment is recommended.
- Use proper local, regional or general anaesthesia.
- Maintain sterile field throughout the procedure and minimize salivary contamination to the implants and surgical site.
- Thoroughly prepare the surgical site and debride the defect of any granulomatous tissue.
- Completely cover the defect area with the membrane.
- Make every attempt to obtain tension-free primary closure over the material.

Preparation of the membrane

- Start the preparation of the membrane when starting the surgical procedure.
- The membrane is supplied in one standard size (30 x 40 mm) together with three different shaping templates (see Figure 1).
- The membrane and templates are packed in a plastic blister. After opening the blister, remove the templates from the packaging and pour all the plasticizer from bottle onto the membrane (cavity A) (see Figure 2). The membrane must be fully immersed in the plasticizer for 20-30 seconds (longer immersion time may reduce mechanical strength).
- Pick-up the membrane using forceps and hold it vertically to allow the excess plasticizer to drain away. Wipe any remaining drops of plasticizer from the lower edge of the membrane. Thereafter, place the membrane in cavity B of the blister (see Figure 3).
- Allow the plasticizer to diffuse in the membrane for at least 10 minutes. The membrane becomes malleable (see Figure 4).
- Trim the membrane to the desired shape using sterile surgical scissors. The preshaped templates are provided as an aid to shaping. Discard the templates after use. DO NOT implant the templates! Take care to round any sharp edges or corners of the trimmed membrane before implantation to prevent irritation of the soft tissue over the membrane.
- At normal room temperature and humidity, the membrane stays malleable for several hours. When exposed to water, as in mouth, the membrane starts to stiffen and will reach maximum stiffness after 30 minutes. The stiffness lasts as long as the membrane remains intact. The stiffening of the membrane can be accelerated by copious flushing with sterile water prior to implantation or once the membrane has been secured into position.

Implantation and fastening of the membrane

- The membrane should always be stabilized, i.e., by securing with tacks or sutures.

- The tacks are contained in a plastic dispenser. Open the dispenser carefully and place it securely on a firm flat surface before fastening the tacks to the tack applicator to avoid dislodging any tacks. Take the tack applicator and press it at an angle of 90 degrees, thus fastening the tack to either the straight or angled end of the applicator.
- Place the shaped membrane in the desired position.
- Drill a hole through the membrane into the bone.
- While keeping the membrane immobile, take the loaded tack applicator and press the tack firmly into the hole. When using the angled end apply finger pressure to the head of the applicator and not to the shaft.
- Turn the tack applicator to one side to release it from the tack.
- Also commercially available resorbable sutures with a non-cutting needle can be used to fasten the membrane.

Postoperative reminders

- As with any surgical procedure, careful postoperative management and patient monitoring is important for optimal healing.
- Provide the patient with detailed instructions for postoperative care.
- Good oral hygiene maintenance is important. The patient should be instructed to rinse with an anti-microbial mouth rinse (e.g. chlorhexidine) between meals for six weeks postoperatively. The patient should not brush or floss teeth in and around the treated area during this time or until the clinician feels that sufficient healing has occurred.
- Antibiotic therapy at the discretion of the clinician.
- Exposed material should be frequently monitored and may be removed postoperatively any time at the discretion of the clinician using sterile surgical scissors.

Evaluation of results

Sites treated with the **INION GTR™ BIODEGRADABLE MEMBRANE** should not be probed for at least 6 months. Gain in attachment level, decreased probing pocket depth, and overall health of the site are effective measurements of determining success of the procedure. Radiographs can be taken to evaluate bone fill postoperatively.

WARNINGS

- Incorrect selection, placement, positioning, and fixation of the implant can cause subsequent undesirable results. The surgeon should be familiar with the devices, the method of application and the surgical procedure prior to performing the surgery.

PRECAUTIONS

- Instruments are available to aid accurate implantation of the **INION GTR™ BIODEGRADABLE MEMBRANE**. Surgical instruments are subject to wear with normal usage and may break. Surgical instruments are only to be used for their intended purpose. All instruments are to be regularly inspected for wear and damage. Use only the **INION GTR™** instruments.
- The patients should have established their willingness and ability to perform adequate oral hygiene prescribed by the clinician.
- Smoking may affect outcomes following periodontal surgery. Treatment of patients who smoke is at the discretion of the clinician.
- Repeated or prolonged contact with NMP may cause skin irritation. If skin contact does occur, the affected area should be washed thoroughly with soap and water for at least 15 minutes.
- As NMP is an eye irritant, if it enters the eye, flush with water for at least 15 minutes.
- Overexposure to NMP solvent vapours can bring about nausea, headache or dizziness.
- DO NOT use for unintended applications! Proper function (i.e. effectiveness and safety) of these implants can not be guaranteed in case of off-label use.

SPECIAL PATIENT POPULATIONS

The **INION GTR™ BIODEGRADABLE SYSTEM** has not been clinically tested in pregnant women.

ADVERSE EFFECTS

Possible complications with any periodontal surgery include thermal sensitivity, gingival recession, flap sloughing, resorption or ankylosis of the treated root, some loss of crestal bone height, perforation or abscess formation, infection, pain, swelling, inflammation, gingival irregularities, and complications associated with the use of anaesthesia. Implantation of foreign materials can result in an inflammatory response or allergic reaction. Transient local fluid accumulation may occur in sterile circumstances.

STERILITY

The **INION GTR™ BIODEGRADABLE SYSTEM** implants have been sterilized with ionizing irradiation. Use only devices that are contained in unopened and undamaged packages. For single use only. Do not re-sterilize. If the device is re-sterilized, the properties of the implant will be affected, and the safety and efficacy cannot be guaranteed. Do not use implant beyond the expiration date on the label.

STORAGE

Store at room temperature (15 to 30°C / 59 to 86°F) at a normal relative humidity. Product should not exceed a maximum temperature defined on the indicator. Do not use if the indicator dot is black.


CAUTION

Federal law (USA) restricts this device to sale by or on the order of a licensed physician or dentist.

MANUFACTURER

INION Oy
Lääkärintie 2
33520 Tampere
Finland
Tel. +358 10 830 6600
Fax +358 10 830 6601
info@inion.com
www.inion.com

Inion® is a registered trademark of Inion Oy.

1

2

3

4

SYMBOLS USED IN LABELING

MD	Medical device	STERILE R	Sterilized using irradiation
	Manufacturer	LOT	Lot number
	Date of manufacture	REF	Catalogue number
	Keep dry		Caution
	Do not re-sterilize		Do not re-use
	Use-by date		Storage temperature limits
	Double sterile barrier system		Do not use if package is damaged and consult instructions for use
UDI	Unique device identifier		Consult instructions for use

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Inion® é uma marca comercial registada da Inion Oy.